



BLOOD STORAGE for EMS



TEMPERATURE

Maintains blood temperature of 1°C-6°C, meeting AABB standards



PORTABILITY & ACCESSIBILITY

Easily transported by vehicle and carried by hand

Alerts and real-time temperature logging through the mobile app



POWER

110V / 12V Vehicle power

Internal charger charges battery while maintaining temperature

24h+ Battery Life



SAFETY

Safety Certified

Helicopter Air Ambulance Qualified



FDA Listed
Blood Refrigerator

