


LifeFlow^{PLUS}
Blood & Fluid Infuser

QinFL^W
NO MORE COLD BLOOD & FLUIDS



A New Paradigm for Rapid Transfusion

Used together, LifeFlow PLUS and QinFlow Warrior lite offer a transformative alternative to traditional rapid infusers for fast, controlled delivery of warmed blood or crystalloids.



Portable, available
when and where
you need it



Easy to set up and
use



A bridge to your
rapid infuser



Use in ER, OR, ICU,
during transport, and
at point of injury

Product Highlights and Features



Contains no aluminum in the warming pathway

Capable of delivering 500ml whole blood at 38°C in less than 3 minutes

Straightforward training and set-up

Works with smaller-gauge catheters



Manufacturer	410 Medical	QinFlow
		
Product Name	LifeFlow PLUS Blood & Fluid Infuser	Warrior lite Blood and Fluid Warmer
Product SKU #	FG5005	Q1310S0000
Contact	info@410medical.com	info@qinflow.com

*Infusion speeds are dependent upon specific equipment and user inputs.
Warrior lite (QinFlow) and LifeFlow PLUS (410 Medical) should each be used in accordance to their respective IFU's.



Scan for Demonstration Video