



MERET oxygen regulators incorporate the most advanced technology to deliver the highest level of safety, reliability, durability, and precision flow.



RESPOND TO THE CALL

949.502.3210
info@meretusa.com
MERETUSA.COM




**OXYGEN
REGULATORS**

DISTRIBUTED BY 

Life-Assist 
Helping Heroes Save Lives

MAGNUM SERIES ALL-BRASS LIFETIME WARRANTY REGULATORS



BRUUT SERIES ALL-BRASS REGULATORS



MERET CGA REGULATOR OPTIONS

REGULATOR	MODEL	FLOWS (L/MIN)	OUTLETS	CONSTRUCTION	WARRANTY
MAGNUM	EMSREG8725	0-25	2 PTO + Barb	Solid 360 Brass, Chrome Plated	LIFETIME
BRUUT	EMSREG8725-xxx	0-25	0-2 PTO + Barb	Solid 360 Brass, Chrome Plated	2 yr
STANDARD	AREG8725-B2D-BU	0-25	2 PTO + Barb	Aluminum/ Brass Internal, Anodized	2 yr
MINI	AREG8715ME-BU	0-15	Barb Only	Aluminum/ Brass Internal, Anodized	2 yr
PRESSURE REDUCER	PRREG54-2D	50 PSI Preset	2 PTO	Solid 360 Brass, Chrome Plated	2 yr

All regulators meet or exceed standards ASTM G175, CGA 870 V-1, CGA G-4.1-2009, ISO 15024-1, EN ISO 13485:2016

BRUUT



- All brass unibody, metal flow knob
- 0-2 PTO and barb

STANDARD



- Aluminum outer body, brass internal
- Multiple flow options +2 PTO

MINI



- Smallest O₂ reg at only 4.0" in length, 0.7 LB

PRESSURE REDUCER



- All brass, compact
- 2 PTO w/ check valves

BRASS VITON WASHER SEALS QTY 50



- FDA recommended for proper seal between a high pressure O₂ cylinder/regulator, offering the safest the option for sealing surface
- Part Number AREG-100WB50