

WATER-JEL BURN CARE THE CLEAR CHOICE

“We keep using Water-Jel because it makes our patients feel better”

— Keith Wesley, M.D., Medical Director, Health East Medical Transportation, St. Paul, MN

“Just used Water-Jel on a welder with a 2nd degree burn. Worked beautifully, as always!”

— Lee Ann Santee, RN, COHN/CM | Project Nurse July 18, 2017



NEW Clear Water-Jel cools, soothes, covers and protects burn injuries.

On the job, at the point of injury, in EMS and in emergency departments, Water-Jel helps providers do more for burns.

- ◆ Provides a consistent Standard of Care for burns
- ◆ Delivers cooling comfort of water – in a pre-saturated non-woven pad
- ◆ Soothes pain on contact
- ◆ Compatible with burn care protocols – especially for more common burns of <15% BSA
- ◆ Specially designed dressings for facial burns and hand burns available
- ◆ Utilizes a clear, smooth, lightweight, water-based gel
- ◆ Allows for visual assessment of burn injuries at any time
- ◆ No pre-wetting or re-wetting required
- ◆ Odor-free.
- ◆ Contains hyaluronic acid and MikroKill
- ◆ Five Year shelf life

Water-Jel burn dressings help absorb excess heat out of the burn site to quickly soothe burn injuries and relieve pain. They cool by convection not evaporation to deliver controlled cooling that is safe and effective for any depth and extent of thermal, electrical friction and radiation burns.

Each burn dressing manages an injury size of 10% or less BSA (shown on the package) so they meet most national, regional and local EMS guidelines. Available in six different sizes, including special Hand and Facial burn dressings, their long shelf life makes them perfect for emergency department and mass casualty incident preparedness and use. Fire Blankets come in three sizes, in orange canisters and foil pouches, to extinguish flames and smoldering clothes while providing comfort. Note: two larger size Fire Blankets do not contain hyaluronic acid.

Water-Jel makes a real difference for burn patients.



Patient with Partial-thickness Steam Burn



Application of Water-Jel Facial Dressing



40 minutes later



Hand Burn Dressing

Take care of serious hand burns – they are very painful and often difficult to manage. Specially designed tabs weave between fingers and wrap the thumb and hand to deliver cooling to all hand and finger surfaces.



Facial Burn Dressing

Facial Burns are critical burn injuries, very painful and difficult to manage. Specially designed dressing features pre-cut holes for the eyes, nose and mouth to deliver cooling to the entire face, forehead, ears, chin and neck.



Add these dressings to your trauma kits or choose from these burn kits.



Item#	Description	Critical Burn Kit for Hands & Face BK-CK-HA bag size:13" x 8" x 1.75"	Emergency Burn Kit Small EBK1-HA bag size:13"x 8"x 1.75"	Emergency Burn Kit Large EBK2-HA bag size:15.4" x 12.5" x 3.5"
B0206-60	2"x6" Burn Dressing			
B0404-60	4"x4" Burn Dressing	2	2	3
B0416-28	4"x16" Burn Dressing		2	1
B0818-60	8"x18" Burn Dressing			1
B1216-20	12"x16" Face Burn Dressing	1	1	1
B0820-20	8"x20" Hand Burn Dressing	2		
G3630P-4	36"x30" Burn Wrap			1

Plastic Burn Kits (Industrial/Welding and Food Service Manufacturing) are also available.
All burn kits contain instruction card, 3" gauze roll and scissors.

For more information, please visit www.Waterjel.com



Water-Jel Technologies L.L.C.
50 Broad Street, Carlstadt, NJ 07072 USA
201-507-8300 waterjel.com waterjelEMS.com

