

NO PPE
REQUIRED!*

CaviWipes™ Bleach

Surface Disinfectant

With our patented chemistry, these multi-purpose cleaning wipes are gentle on surfaces, yet tough on pathogens. They kill *C. diff* spores in 3 minutes and do not require PPE*.

They also meet the CDC recommended 1:10 bleach dilution and are compatible with most medical device equipment**. So when it comes to disinfecting hard non-porous surfaces, leave the gloves behind and reach for CaviWipes Bleach.



3-Min Universal Contact Time against 43 relevant pathogens including *C. Diff*, TB, *C. auris*, Norovirus, and 8 MDROs (Multidrug Resistant Organisms).



Fully qualified for EPA Emerging Viral Pathogen (EVP Level III) Claim for all virus types (enveloped, large and small non-enveloped viruses) to meet current and future infection prevention needs.



No Signal Word - No PPE required*.



Minimal Residue, good cleaner and high substrate quality.



Mild odor.



Compatible with most medical equipment.

DISTRIBUTED BY 

Life-Assist
Helping Heroes Save Lives 

Metrex™
INFECTIOUS PREVENTION

CaviWipes Bleach demonstrates effectiveness against the following 43 relevant microorganisms in 3 minutes:

Bacteria

Clostridium difficile spores
Staphylococcus aureus
Pseudomonas aeruginosa
Salmonella enterica
Burkholderia cepacia
Escherichia coli
Klebsiella pneumoniae
Legionella pneumophila
Listeria monocytogenes
Serratia marcescens
Streptococcus pyogenes
Neisseria gonorrhoeae

Drug Resistant Bacteria

Enterobacter cloacae
Multidrug Resistant *Acinetobacter baumannii* (MRAB)
Carbapenem-resistant *Klebsiella pneumoniae* (CRKP) (CRE)
Methicillin Resistant *Staphylococcus aureus* (MRSA)
Penicillin Resistant *Streptococcus pneumoniae*
Vancomycin Resistant *Staphylococcus aureus* (VRSA)
Vancomycin Resistant *Enterococcus faecalis* (VRE)
ESBL *Escherichia coli*
Mycobacterium tuberculosis var: *bovis* (BCG) (TB)

Mold, Mildew, Fungi

Candida auris
Candida albicans
Aspergillus brasiliensis
Trichophyton interdigitale

Viruses (Non-Enveloped)

Adenovirus
Hepatitis A virus
Norovirus (as Feline Calicivirus)
Poliovirus
Rhinovirus
Enterovirus (EV-D68)
Rotavirus

Viruses (Enveloped)

Herpes simplex virus type 1
Herpes simplex virus type 2
Human Coronavirus
Influenza A virus Strain A/Hong Kong 8/68 (H3N2)
Influenza B virus Strain B/Hong Kong/5/72
Hepatitis B virus (HBV)
Hepatitis C virus (HCV)
Human Immunodeficiency virus (HIV)
Respiratory syncytial virus (RSV)
Measles virus
SARS-CoV-2



Ordering Information

Part No.	Description
13-9100	CaviWipes Bleach (6" x 10.5") – 90 towelettes per canister
13-1175	CaviWipes Bracket

*Special cleaning and PPE instructions Prior to Disinfection against *Clostridium difficile* spores, HIV-1, HBV, HCV.
**See medical device equipment manual for compatibility information.

To learn more about CaviWipes Bleach scan the QR code or visit metrex.com

