

FORA® *Test N'GO Advance Pro*

MULTI-FUNCTIONAL MONITORING SYSTEM (FOR MULTIPLE-PATIENT USE)

FEATURES

- Gold advance strip technology
- Tests blood glucose and β -Ketone
- Bluetooth connectivity
- LCD backlight & Strip indication light

SPECIFICATIONS

Blood Glucose Test Volume: 0.5 μ L
Blood Glucose Reaction Time: 5 seconds
Measurement Range: 10 ~ 600 mg/dL
Memory Capacity: Stores 1,000 test results
Blood Ketone Test Volume: 0.8 μ L
Blood Ketone Reaction Time: 10 seconds
Measurement Range: 0.1 ~ 8.0 mmol/L
Power Source: 1 x AAA



Multi-Functional Blood Monitors for Health Care Professionals

The **FORA Test N'GO Advance Pro** multi-functional blood monitors are used in professional settings such as EMS, long-term care facilities or hospitals to measure blood glucose and ketone levels. Several features to ensure a proper care for patients at bedside include strip ejector to avoid unnecessary contact, GDH enzyme for venous, or capillary blood, and Bluetooth connectivity to apps and EHR to share with other physicians. Using the most advanced Gold Strips Technology, our devices guarantee measurement results that meet the highest FDA standards. Individually foil wrapped packages for the strips are also available for one-time use to prevent contamination.

DISTRIBUTED BY 
Life-Assist
Helping Heroes Save Lives 