

**NEW**

# microdot<sup>®</sup> PRO

## Blood Glucose Monitoring System

For use by Health Care Professionals



The microdot<sup>®</sup> PRO Blood Glucose Monitoring System incorporates advanced algorithm and sensing technologies to dramatically improve the testing result on blood samples with extreme hematocrit range to ensure more accuracy in each test



**Small Blood Sample**  
Sample size .8  $\mu$ l



**High Visibility Meter**  
and protective sleeve



**Strip Ejection**



**GDH-FAD**  
for less interference



**Fast Results**  
8 second test



**Autocode**



**Venous & Capillary**  
whole blood approved

**microdot<sup>®</sup> PRO kit includes:** microdot<sup>®</sup> PRO Blood Glucose Meter, High Visibility Protective Sleeve, and Large or Small Pouch

### Specifications

**Measurement Range :** 20-600 mg/dL (1.1 ~ 33.3 mmol/L)

**Testing Time (seconds) :** 8 seconds

**Coding Method :** Automati

**Hematocrit Range (%) :** 10-70%*c*

**Sample Volume ( $\mu$ L) :** >0.8  $\mu$ L

**Operating Temperature ( $^{\circ}$ C) :** 5 $^{\circ}$ C ~ 45 $^{\circ}$ C

**Memory Capability :** 1000 test results

Ask for *microdot*<sup>®</sup>

DISTRIBUTED BY 

**Life-Assist**

Helping Heroes Save Lives 

Call or visit our website for our full line of products [microdotcs.com](http://microdotcs.com)

toll free: 877.374.4062 | [info@microdotcs.com](mailto:info@microdotcs.com)

**CS** CAMBRIDGE SENSORS USA, LLC

**NEW**

# microdot<sup>®</sup> PRO

## Blood Glucose Monitoring System

The microdot<sup>®</sup> PRO Professional Blood Glucose Monitoring System includes:  
microdot<sup>®</sup> PRO Blood Glucose Meter, Test Strips, Control Solutions, Protective Sleeve and Pouch.



### No INTERFERENCE from Over 25 Substances at Therapeutic Levels

- |                 |                   |                 |              |                |
|-----------------|-------------------|-----------------|--------------|----------------|
| • Acetaminophen | • Acetylsalicylic | • Maltotetraose | • Icodextrin | • Cholesterol  |
| • L-DOPA        | • Methyldopa      | • Xylitol       | • Fructose   | • Bilirubin    |
| • Tolbutamide   | • Tetracycline    | • Xylose        | • Hemoglobin | • Triglyceride |
| • Dopamine      | • Heparin         | • Gentic acid   | • Creatinine | • Galactose    |
| • Ibuprofen     | • Maltotriose     | • Pralidoxime   | • Uric Acid  | • Glutathione  |

Ask for *microdot*<sup>®</sup>

microdot<sup>®</sup> PRO Glucometer: Product Number: 400-01

microdot<sup>®</sup> PRO Blood Glucose Test Strips: 50 Count: 400-50 | 25 Count: 400-25

microdot<sup>®</sup> PRO Control Solutions: 400-02

Large Pouch: Product Number: 401-LG | Small Pouch: Product Number: 401-SM

Call or visit our website for our full line of products [microdotcs.com](http://microdotcs.com)

toll free: 877.374.4062 | [info@microdotcs.com](mailto:info@microdotcs.com)



**CS** CAMBRIDGE SENSORS USA, LLC