



FOLDABLE PATIENT TRANSPORT CHAIR WITH TRACKS

4BELL STAIR+

FEATURES

- **Frame in aluminium** makes the chair lightweight and resistant at the same time
- **Two reinforced tracks** help going down stairs in a safe and controlled way
- **An extendable handhold** on the back offer greater control when sliding the chair down stairs
- **Two telescopic handles at the front and two foldable ones behind the backrest** help operators to lift the chair in case of obstacles
- **Rear wheels' locking mechanism** makes sure the chair does not slide on a flat surface
- **Semirigid seat and soft backrest** are built to increase patient comfort
- **Foldable footrest** increases patient stability
- **Easy to clean PVC seat and backrest** are resistant to chemicals
- **Pivoting front wheels** makes the chair more maneuverable.



| | |
|--------------|---------------------|
| Model | 4BELL STAIR+ |
| Frame colour | |
| Seat colour | |
| Code | ST10501 |

4BELL STAIR+ is an EMS foldable transport chair designed to transfer or evacuate patients with mobility limitations from multi-storey buildings. Its two tracks make the chair suitable for a controlled stairway descent. Lifting the chair may not be required.

It is built in accordance to the specifications set by UNI EN 1865-4 and **is compliant with UNI EN 1789 when anchored to its dedicated fixation system 4BELL Stair+ Max.**

The frame of the chair is made in aluminium for increased resistance. Its total weight is 14,7 kg (32.4 lb).

180 kg - 397 lb



Load Capacity



Stair Tracks



Multi-height Extension Handle



Foldable Chair

TECHNICAL DATA

| | |
|---|---------------------------------------|
| Length with handles retracted / extracted | 76 / 112 cm - 29.9 / 44 in |
| Length with tracks and telescopic handles extracted | 128 cm - 50.3 in |
| Width | 49 cm - 19.2 in |
| Height | 94 cm - 37 in |
| Height with extendable handle extracted | 147,5 cm - 58.07 in |
| Dimensions when folded | 49 x 94 x 29 cm - 19.3 x 37 x 11.4 in |
| Weight | 14,7 kg - 32.4 lb |
| Load capacity | 180 kg - 397 lb |
| Front wheels | ø 10 cm - 3.93 in |
| Rear wheels | ø 15 cm - 5.9 in |
| Materials | Al, Steel, PU, PVC |

THE PRODUCT INCLUDES

Restraint belts.

ACCESSORIES



ST10455B

4BELL STAIR+ MAX
 10G tested fastening system

- It is installed in ambulances, special vehicles or airplanes
- It helps securing the 4BELL STAIR+ in foldable position to either the wall or the floor.



| | |
|---------|------------------------|
| ST10452 | Seat XL |
| ST10450 | Armrests |
| ST10451 | Rear auxiliary handles |
| ST10454 | Headrest 4Bell STAIR+ |

Rev.0 (14/03/2022)

©Copyright - This document is property of SPENCER ITALIA S.r.l and can't be reproduced, used and disclosed without the written approval of SPENCER ITALIA SRL.

