

The Meret logo is displayed in a bold, black, sans-serif font. The letter 'M' is stylized with a horizontal bar extending to the left and a vertical bar extending downwards. The word 'MERET' is positioned to the right of the 'M'. A registered trademark symbol (®) is located at the end of the word. The logo is set against a bright yellow rectangular background.

MERET®

A black tactical bag, the V.E.R.S.A.™ PRO X, is the central focus. It features a MOLLE webbing pattern on the front and side pockets, a top handle with two straps, and a main zippered compartment. The bag is shown in a well-lit industrial or warehouse setting. In the background, a large American flag is mounted on a wall, and the structural beams of a ceiling are visible. The bag is positioned in the lower half of the frame, with the text 'V.E.R.S.A.™ PRO X' overlaid on a black horizontal band across its middle.

V.E.R.S.A.™ PRO X

V.E.R.S.A.™ PRO X

Part #: M8127TB
Dimensions: 18.5" x 13" x 10.5"
Weight: 8lbs 1oz



PROFESSIONAL USE

The most versatile backpack for Advanced Life Support (ALS) and Basic Life Support (BLS) applications

FEATURES

1. Ergonomic heavy-duty top carry handle
2. Large EZ-Pull pullers allow easy access
3. Large front tape plate area for ID customization
4. Interior and Exterior Infection Control Material allows for quick decontamination of bloodborne pathogens and other potentially transmissible substances
5. Hideaway zippers secures to all MERET module zipper ready bags
6. Black matte finish metal hardware
7. Double stitching on all stress points
8. Laser cut MOLLE panel for **OPTIONAL** modules or accessories

OPTIONS

9. Shown with **OPTIONAL** XTRA FILL™ PRO X module
10. Ballistic Armored Panel



Modules and supplies shown in the images above are not included.

